
MICROPUMP SYSTEM

CORRECTIONS TO THE USER'S MANUAL

⚠ Important Safety Information / Wichtiger Sicherheitshinweis / Informations importantes sur la sécurité / Informazioni di sicurezza importanti / Información importante de seguridad / Informações importantes de segurança / Informações de segurança importantes / Belangrijke veiligheidsinformatie / Viktige sikkerhedsoplysninger / Viktig informasjon om sikkerhet / Viktig säkerhetsinformation / Tärkeitä turvallisuustietoja / Ważne informacje dotyczące bezpieczeństwa / Důležité bezpečnostní informace / معلومات هامة للسلامة

Chapter 4.1 Recommended Infusion Sites

The note correctly reads as follows:

NOTE

- Do not attach the pump holder to body sites with scars, birthmarks, moles, tattoos, injuries, bruises, or rashes.
- Choose a flat infusion site that is not too curved so that the adhesive pad of the pump holder adheres well.
- Choose an infusion site that is easy to reach and see.
- A new infusion site must be at least 5 cm away from the last infusion site.
- Before you attach the pump holder to the body, the infusion site must be completely dry.
- Do not touch the introducer needle or the cannula of the cannula assembly.
- If the pump holder frequently detaches from the skin, consult your healthcare professional to find a method that may improve adhesion to the skin.
- If the infusion site becomes inflamed or localised skin reactions occur (for example, allergic reaction or eczema), then replace the infusion assembly immediately and choose a new infusion site.

Chapter 4.2.1 Attaching the Infusion Assembly to the Body

The new warning and steps 9 and 10 correctly read as follows:

WARNING

Risk of hyperglycaemia (high glucose level)

Improper handling of the insertion device may lead to leakage. This may lead to hyperglycaemia.

Carefully follow the instructions for using the insertion device as described in the steps below.

- 9 Holding the skin taut, press the insertion device firmly against the selected infusion site without tilting, twisting or wiggling. Apply force only in the direction of the infusion site as shown in the illustration.
- 10 While holding the insertion device firmly in place, press the blue release button to insert the cannula under the skin.

Chapter 4.2.6 Attaching the Micropump

The warning, note and step correctly read as follows:

WARNING

Risk of hyperglycaemia (high glucose level)

- Deformations, cracks, or frequent and incorrect insertion of the micropump into the infusion assembly may cause the micropump system to become leaky. Before attaching the micropump to the infusion assembly, check the micropump and the pump holder for damage.
- Take special care not to bend the reservoir needle while attaching the micropump to the infusion assembly. A bent reservoir needle can cause the micropump system to become leaky.

NOTE

Check the site on your body with the adhesive pad of the infusion assembly at least once a day.

- 2 Keep the micropump straight and parallel to the pump holder. Turn the micropump carefully in place. You may apply light pressure on the pump shield to make sure the reservoir needle stays inside the cannula head of the pump holder. Try to avoid tilting movements.
After one eighth turn (approx. 45 degrees), the micropump clicks into the front and rear hooks on the pump holder.

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