en	Instructions for use/Technical description Silver Graft straight tubes, bifurcations / Silver Graft Helix straight tubes and axillo-bifemorals
USA	Note for U.S. users This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusaifus.com . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.
de	Gebrauchsanweisung/Technische Beschreibung Silver Graft gerade Röhrchen, Bifurkationen / Silver Graft Helix gerade Röhrchen und Axillo-Bifemorale
fr	Mode d'emploi/Description technique Silver Graft tubes droits, bifurqués / Silver Graft Helix tubes droits et axillo-bifémoraux
es	Instrucciones de manejo/Descripción técnica Tubos rectos, bifurcaciones Silver Graft / tubos rectos y axilobifemorales Silver Graft Helix
it	Istruzioni per l'uso/Descrizione tecnica Tubi rettilinei Silver Graft, biforcazioni/tubi rettilinei Silver Graft Helix e axillo-bifemorali
pt	Instruções de utilização/Descrição técnica Tubos retos, bifurcações Silver Graft / Tubos retos e axilobifemorais Silver Graft Helix
nl	Gebruiksaanwijzing/Technische beschrijving Silver Graft rechte grafts, bifurcaties / Silver Graft Helix rechte grafts en axillobifemorals
da	Brugsanvisning/Teknisk beskrivelse Silver Graft lige rør, bifurkationer/ Silver Graft Helix lige rør og axillo-bifemorale
no	Bruksanvisning/Teknisk beskrivelse Silver Graft rette rør, bifurkasjoner / Silver Graft Helix rette rør og axillo-bifemorale
sv	Bruksanvisning/Teknisk beskrivning Silver Graft raka rör, förgrenade/Silver Graft Helix raka rör och axillo-bifemorala
fi	Käyttöohjeet/Tekninen kuvaus Silver Graft suorat putket, bifurkaatiot / Silver Graft Helix suorat putket ja aksillo-bifemoraalit
et	Kasutusjuhend/Tehniline kirjeldus Silver Grafti sirged voolikud, hargmikud / Silver Graft Helixi sirged voolikud ja aksillobifemoraalid
lv	Lietošanas instrukcija/Tehniskais apraksts Silver Graft taisnas caurules, divvirzienu / Silver Graft Helix taisnas caurules un aksilu bifemorālas
lt	Naudojimo instrukcija/Techninis aprašymas Silver Graft tiesūs vamzdeliai, bifurkacijos / Silver Graft Helix tiesūs vamzdeliai ir aksillobifemoraliniai šuntai
ru	Инструкция по применению/Техническое описание Прямые трубки, бифуркационные протезы Silver Graft / прямые трубки и подмышечно-бедренные протезы Silver Graft Helix
cs	Návod k použití/Technický popis Rovné trubičky, bifurkace Silver Graft / rovné trubičky a axilo-bifemorální protézý Silver Graft Helix
pl	Instrukcja użytkowania/Opis techniczny Rurki proste Silver Graft, dwudzielne / rurki proste Silver Graft Helix i pachwinowo-dwuudowe
sk	Návod na použitie/Technický opis Silver Graft Rovné rúrky, rozdvójky / Silver Graft Helix rovné rúrky a axilo-bifemorálky
hu	Használati utasítás/Műszaki leírás Silver Graft egyenes csövek, bifurkációk / Silver Graft Helix egyenes csövek és axillo-bifemorálisok
sl	Navodila za uporabo/Tehnični opis Ravne cevke, dvokrake Silver Graft/ravne cevke in aksilo-bifemoralne cevke Silver Graft Helix
hr	Upute za upotrebu/Tehnički opis Silver Graft ravne cijevi, račvaste cijevi / Silver Graft Helix ravne cijevi i aksilarno-bifemoralne
ro	Manual de utilizare/Descriere tehnică Tuburi drepte, bifurcații/Silver Graft / tuburi drepte și axilo-bifemorale Silver Graft Helix
bg	Инструкции за употреба/техническо описание Прави тръби и бифуркации Silver Graft / Прави тръби и аксило-бифеморални протези Silver Graft Helix
tr	Kullanım Kılavuzu/Teknik Açıklama Silver Graft düz tüpler, bifurkasyonlar / Silver Graft Helix düz tüpler ve aksillo-bifemoraller
el	Οδηγίες χρήσης/Τεχνική περιγραφή Silver Graft ευθείς σωλήνες, διακλαδώσεις / Silver Graft Helix ευθείς σωλήνες και μασχαλοδιμηνιαίοι
br	Instruções de uso/Descrição técnica Tubos retos, bifurcações Silver Graft / tubos retos e axilobifemorais Silver Graft Helix

AESCULAP® Silver Graft

Silver Graft straight tubes, bifurcations / Silver Graft Helix straight tubes and axillo-bifemorals

1 About this document

Note

General risk factors associated with surgical procedures are not described in these instructions for use.

Note

Instructions for use and further information about B. Braun / AESCULAP products can be found on the B. Braun eFU website at eifu.bbraun.com.

The term Silver Graft refers to all variants of this product class (i.e. straight tubes, bifurcations / Silver Graft Helix straight tubes and axillo-bifemorals) if not explicitly stated otherwise.

2 Clinical use

Note

The summary of safety and clinical performance (SSCP) of the product is available in the European database on medical devices (EUDAMED) at ec.europa.eu/tools/eudamed

The SSCP can be found using the following Basic UDI-DI:

- 40392390000024322E
- 40392390000024332G

2.1 Product description

Silver Graft prostheses are warp-knitted, double-velour vascular prostheses made of polyethylene terephthalate (polyester) for reconstructive procedures on the aorta and arteries in the peripheral area. They are impregnated with absorbable modified bovine gelatin for temporary sealing. The outer surfaces of the prostheses are coated by a metallic silver layer to inhibit microbial adhesion. Temporary sealing and inhibition of microbial adhesion are confirmed by product tests as well as preclinical in vivo and in vitro tests.

Silver Graft prostheses are available in different lengths and diameters.

Additional information is available from B. Braun at <https://www.bbraun.com> or from your local B. Braun / Aesculap agency.

2.1.1 Mode of action

The functional properties of Silver Graft are based on the principle that a highly porous textile construction is temporarily sealed by a biologically degradable mass, which does not possess any independent pharmacological properties. This modified bovine gelatin impregnation reduces the initial porosity of the textile prosthesis at the time of implantation of approximately 1 500 ml/min/cm² (Silver Graft Helix) or 1 900 ml/min/cm² (Silver Graft) at a pressure of 120 mm Hg to 0 ml/min/cm². From outside the impregnation of the prosthesis is enzymatically degraded and absorbed within 4–6 weeks after implantation, as shown in a publication by Ueberueck et al. in 2005. During this time endogenous tissue incorporates the prosthesis into the surrounding tissue and takes over the sealing function. On the inner side of the prosthesis a cell layer is formed, referred to as pseudo neointima.

Silver Graft Helix is wrapped with an external polypropylene thread around the outside of the prosthetic wall to protect the prosthesis against mechanical pressure from outside and increase the kink-stability.

The silver ions on the graft surface have an antimicrobial (bacteriostatic and bactericidal) effect. Therefore, Silver Graft has a significant reduction of microbial adhesion to the surface of the textile prosthesis confirmed by preclinical in vivo and in vitro tests. Hereof a slow release of the silver ions ensures long-time antimicrobial properties. Investigations in cell cultures as well as in animals show the biocompatibility of Silver Graft.

2.1.2 Product range

The product-range consists of straight tubes, straight tubes helix supported and bifurcations of various lengths (15 cm to 90 cm) and diameters (6 mm to 24 mm) as well as axillo-bifemoral prostheses.

All forms and sizes of Silver Graft vascular prostheses are supplied in a double-peel-packaging and are closed in a cardboard box. Box content is 1 piece.

2.1.3 Materials

Silver Graft components are:

- Polyethylene terephthalate (polyester)
- Absorbable modified bovine gelatin, concentration of 10 to 25 wt% for Silver Graft and 8 to 17 wt% for Silver Graft Helix
- Glycerol
- Hexamethylene diisocyanate (crosslinker for gelatin)
- Metallic silver, concentration of 0.17 to 0.40 wt% for Silver Graft and 0.16 to 0.35 wt% for Silver Graft Helix
- Polypropylene
- Palladium
- Titanium

2.2 Areas of use and limitations of use

2.2.1 Intended use

Reconstructive procedures on the aorta and in the peripheral area requiring a prosthesis with at least 6 mm in diameter. Axillofemoral and –bifemoral reconstructions are also indicated.

Silver Graft can be implanted prophylactically as a protection against vascular graft infection. It is designed for long-term protection avoiding adhesion of microbes to the surface of the textile prosthesis.

2.2.2 Indications

Note

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

Silver Graft is intended for use as vascular prosthesis for replacement or bypass of diseased vessels in patients suffering from aneurysmal or occlusive disease. The supported helix variant is particularly intended for use in locations that are exerted to increased mechanical stress. E.g. this is the case for bypasses that are located directly under the skin or for above the knee.

2.2.3 Absolute contraindications

- DO NOT use Silver Graft in patients
 - who have had a previous incidence of heparin-induced-thrombocytopenia type II (HIT II).
 - for any configuration of coronary artery bypass procedures, cerebral reconstruction, blood access fistula (e.g. hemodialysis) or pulmonary shunting
 - with known sensitivity to silver, palladium, polyester or materials of bovine origin or
 - as a patch.
- DO NOT use Silver Graft Helix in Children < 10 kg

2.2.4 Relative contraindications

The clinical use of the product depends on the knowledge and experience of the surgeon whose responsibility is to decide whether the product can successfully be implanted, considering the safety and warning information in the instructions for use.

Pregnancy and lactation

No clinical data is available on the usage of Silver Graft in pregnant and breastfeeding women. No contraindications and risks concerning reproduction are known. The treating surgeon should balance the medical advantages of a silver coated vascular prosthesis against the potential risks.

2.3 Risks, adverse effects and interactions

It is known, in connection with the use of vascular prostheses that the following side effects, but not exclusively, can occur:

- Seroma in the area of the prosthesis
- Mechanical injury or tearing of sutures, prosthesis and/or host vessel
- Occlusion/stenosis of the prosthesis
- Thrombosis
- Bleeding
- Hematoma in the area of the prosthesis
- Anastomotic aneurysms
- Perigraft reactions
- Dilatation of the vascular graft
- Embolisms
- Perigraft seroma formation
- Local argyria
- Allergic reactions to gelatin of bovine origin.

The physician should consider the need for intraoperative and/or postoperative anticoagulation therapy based on the pharmacological requirements and medical history of the patient. Effective anticoagulation and antiplatelet therapy should be maintained at a dosage deemed appropriate by the physician.

Careful and continuous medical follow-up is advised so that prosthesis-related complications can be diagnosed and properly managed to minimize danger to the patient.

2.4 Safety information

2.4.1 Clinical user

General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- Use the product only according to these instructions for use.
- Follow the safety and maintenance instructions.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- Keep the instructions for use accessible for the user.
- Do not use damaged or surgically removed implant components.

Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

Notes on surgical procedures

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

It is the user's responsibility to ensure the correct implant size and implantation.

Aesculap is not responsible for complications caused by:

- Incorrect indication or implant selection
- Incorrect surgical technique
- Incorrect combination of implant components, particularly in combination with components of other manufacturers

- Exceeding the limitations of the treatment method or non-observance of essential medical precautions

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

2.4.2 MRI Safety and Compatibility

Magnetic resonance Imaging (MRI) can be safely done in patients with implanted Silver Graft vascular prostheses under the following conditions:

- Static magnetic field: 1.5 T and 3 T
- Maximum spatial field gradient: 6.6 T/m
- Whole body averaged specific absorption rate (SAR): 2 W/kg

2.4.3 Sterility

The product has been sterilized by irradiation (Silver Graft) or ethylene oxide (Silver Graft Helix) and is supplied in sterile packaging.

The safe sterile provision of the product is only guaranteed if the sterile packaging is undamaged and unopened and the use-by date has not passed.

- Store products in their original packaging in a dust-protected, dry and temperature-controlled room.
- Remove products from their original protective packaging only just prior to application.
- Visually inspect the product packaging to ensure that the sterile barrier system is intact.
- Do not use products from open or damaged sterile packaging.
- Do not use the product after its use-by date.

The product is intended for single-use.

- Do not reuse the product.

The reprocessing of the product affects its functionality. Risk of injury, illness or death due to soiling and/or impaired functionality of the product.

- Do not reprocess the product.

2.4.4 Transport and storage

The following environmental conditions apply to the transport and storage of the product:

Temperature	max. 27 °C
-------------	------------

2.5 Patient education

Within the framework of the patient education, the relevant circumstances needed for consent must be explained to the patient in accordance with their level of understanding, pre-existing knowledge and need for information. This includes:

- Diagnosis, procedure and risk clarification
- Operative procedure
- Advantages and disadvantages of the procedure
- All alternative procedures that can be considered

The patient must be properly informed about the procedure and in particular about the following information:

- The patient must undergo medical check-ups of the implant at regular intervals.
- Pains may occur following the treatment.
- The patient must be informed about the limits to the allowable strain on the implant and be given appropriate behavior guidelines. The risks of transgressing these rules must be explained to the patient.

2.6 Application

2.6.1 Documentation

The user shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant and its dimensions
- Location of intraoperative landmarks

Important information concerning the implanted product and the surgical procedure can be noted down on the implant card. The implant card is enclosed with the product.

Note

A printable version of the implant card, a patient information brochure and a guide on how to fill out the enclosed implant card (if applicable) are available online at www.bbraun.com/vascularimplantpassport

- Fill out the implant card and hand it to the patient together with the patient information brochure after the surgical procedure.

Each package contains additional labels showing the designation, article and lot number and – if applicable – the individual serial number of the product.

- Use these labels for documentation in the patient's file (for the hospital) and the implant card (for the patient).

The patient's consent must be documented.

2.6.2 Application notes

Choose the suitable shape of prosthesis and the length and calibre of prosthesis to suit the indication and the anatomical situation.

It is not necessary to seal with endogenous blood as the prosthesis is blood tight.

The suture materials with atraumatic round-body needles commonly used in vascular surgery are suitable for fixing the vascular prosthesis. The unsupported Silver Graft prosthesis should be sutured "under slight tension", i.e. the most to the point until the pleating is flat.

When using Silver Graft Helix, the spirals can be removed at each end to facilitate the anastomoses.

It is possible to apply antibiotics to the prosthesis (e.g. Rifampicin) for prophylaxis according to recommendations of the German Society of Vascular Surgery, DGG. The prosthesis surface, however, has to be completely immersed in the antibiotic solution.

The first choice to treat a vascular graft infection is the surgical intervention. Infected vascular graft material should be completely removed and the principals in septic surgery should be observed. Autologous in situ re-constructions have the lowest risk of reinfection. Silver Graft should be implanted if autologous material is not available or not possible to use given the patient's risk/benefit status.

3 Disposal

⚠ WARNING

Risk of infection due to contaminated products!

- Adhere to national regulations when disposing of the product, its components and its packaging.

4 Symbols on product and packaging

	Follow the instructions for use
	MR conditional
	Medical device
	Contains a medicinal substance
	Contains biological material of animal origin
	Double sterile barrier system
	Sterilization using irradiation (Silver Graft straight tubes and bifurcations)
	Sterilization using ethylene oxide (Silver Graft Helix straight tubes and axillo-bifemoral)
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Date of manufacture
	Size
	Manufacturer's batch designation
	Manufacturer's article number
	Upper limit of temperature
	Do not use if packaging is damaged
	CE marking according to regulation (EU) 2017/745

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